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000978EX-20141209

Pi-Ku-Plast

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Instructions for use

Pi-Ku-Plast

Befor using the product, please read these instructions for use carefully!

The manufacturer will not accept any liability for damage resulting from non-compliance with these instructions for use.

1. Description of the product

Pi-Ku-Plast is a cold-curing brush resin available in various colors for modelling dental structures or restorations.

2. Indication

- Pi-Ku-Plast is used
 - for modelling of secondary elements for all types of attachments and double crowns
 - for stress-free fixation and connection of crowns for soldered restorations and for wax patterns/models
 - for the fabrication of resin dies
 - for copings on resin dies; with the use of Pi-Ku-Plast separating varnish (REF 540 0018 2) as insulating agent

3. Contraindication

Pi-Ku-Plast is not suitable for direct (intraoral) use.

4. Hazard and safety information

Pi-Ku-Plast Cleaner

Signal word:

Danger

Hazard statements:

- H225 Highly flammable liquid and vapour.
- H319 Causes serious eye irritation.
- EUH066 Repeated exposure may cause skin dryness or cracking.

Precautionary statements:

- Prevention:**
 - P210 Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.
 - P233 Keep container tightly closed.
 - P280 Wear protective gloves/protective clothing/eye protection/face protection.

Response:

P303+P361+P353 IF ON SKIN (or hair): Remove/Take off immediately all cantaminated clothing. Rinse skin with water/shower.
P305+P351+P338 IF IN EYE: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

Storage:

P403+P235 Store in a well-ventilated place. Keep cool.

Pi-Ku-Plast Polymer

Signal word:

Warning

Hazard statements:

- H317 May cause an allergic skin reaction.

Precautionary statements:

- Prevention:**
 - P261 Avoid breathing dust
 - P280 Wear protective gloves/eye protection

Response:

P333+P313 If skin irritation or rash occurs: Get medical advice/attention.
P363 Wash contaminated clothing before reuse.
P302+P52 IF ON SKIN: Wash with plenty of soap and water.

Disposal:

P501 Dispose of contents/container in accordance with local/regional/national/international regulations.

Product contains:

dibenzoyl peroxide, methyl methacrylate

Pi-Ku-Plast Monomer

Signal word:

Danger

Hazard statements:

- H225 Highly flammable liquid and vapour.
- H315 Causes skin irritation.

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H335 May cause respiratory irritation.

Precautionary statements:

Prevention:

- P210 Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.
- P233 Keep container tightly closed.
- P261 Avoid breathing/vapours.
- P280 Wear protective gloves/eye protection.

Response:

P303+P361+P353 IF ON SKIN (or hair): Remove/Take off immediately all contaminated clothing. Rinse skin with water/shower.

Product contains:

methyl methacrylate, ethylene dimethacrylate, 2-[2-H-Benzotriazol-2-yl]-p-cresol

5. Storage and durability

The product must be stored in a cool, well-ventilated room.

6. Processing

6.1 Preparation of the object

The plaster dies must be insulated once with Pi-Ku-Plast separating varnish before modelling the resin copings. Apply Pi-Ku-Plast after a drying time of approx. 5 min. The metal surfaces of primary elements should be polished to a high luster and should not have undercuts in order not to render the re-moval of Pi-Ku-Plast more difficult.

6.2 Preparatory work

Fill polymer, monomer and cleaner into the containers included in the Pi-Ku-Plast assortment. The color is included in the mo-nomer. Put the enclosed round colored marking that corresponds to the color of the monomer onto the monomer container.

Please note:

- Polymer and monomer must not be poured back into the container. (Caution: Cleaner must be removed from the brush before the brush is immersed into the monomer.) It can be stored for few days in well-sealed containers.
- Do not mix the cleaner with the monomer. (Caution: the cleaner affects the accuracy of fit, resistance to fracture, removal behavior and hardening).
- If polymer (powder) gets into the working container of the monomer (liquid), the liquid must be replaced. Otherwise the monomer polymerizes within 2 days - even in a well-sealed working container.
- After the use, the container must be immediately sealed again.
- Do not spill monomer - staining will result!

6.3 Processing

Wet the brush with monomer and then take up the polymer with the tip of the brush. Depending on the monomer quantity on the brush as well as the duration of immersion and the immersion depth in the polymer, the quantity to be taken up can be varied. By keeping the resin on the tip of the brush for two to five seconds, a plastic modelling behavior is achieved so that even cusps and ridges can be easily applied. Accordingly, it is possible to take up the exact quantity of Pi-Ku-Plast required and to control the modelling behavior. Contamination of the monomer with polymer or also with the cleaner (wet brush after cleaning) can affect the properties of the resin (increased bonding to the metal, reduced hardness).

Recommendation: The monomer in the working container should be exchanged every day!

Pi-Ku-Plast can be perfectly milled and does not form any gaps when it is applied onto Pi-Ku-Plast material that has already polymerized. Clean the brush before Pi-Ku-Plast has hardened.

Plast material that has already polymerized. Clean the brush before Pi-Ku-Plast has hardened.

Please note: After cleaning, the brush must be dried well! The cleaner must not get into the monomer!

6.4 Processing time

The processing time varies depending on the room temperature. At 23° C it amounts to four to five minutes.

Further processing

Pi-Ku-Plast can be applied in layers. Pi-Ku-Plast can also be applied to hardened Pi-Ku-Plast material without the for-mation of fine gaps. The monomer color can be changed individually at any time; this has no influence on the quality of the pattern/model. All modelling waxes can be applied to Pi-Ku-Plast. Pi-Ku-Plast can be processed and polished with all commercially available rotary instruments. The product is also perfectly suitable for connecting wax patterns since it generates hardly any setting heat. Pi-Ku-Plast does not generate any expansion forces when it is burned, but it collapses; this way damage to the shape prior to the cristobalite crack is avoided.

6.5 Use on patients (direct use)

Pi-Ku-Plast is not an approved medical product (device) and must not be used in the patient's mouth.

7. Symbols

7.1 Symbols according to DIN EN ISO 15223–1

REF Reference number

 Manufacturer

7.2 Symbols according to GHS regulation

 GHS 02

 GHS 07

7.3 Symbols according to NFPA (USA)

 Hazard diamond

8. Additional information

These instructions for use are based on state-of-the-art methods and equipment and our own experience. The product may only be used for the indication described under item 2. The user himself is responsible for processing the product. Liability for incorrect results shall be excluded since the manufacturer does not have any influence on further processing. Any occurring claims for damages may only be made up to the value of our products.

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